

Priz/TPD52 protein (STJ505030)

STJ505030

GENERAL INFORMATION

Product Type	Proteins
Short Description	Priz/TPD52-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7163
Gene Symbol	TPD52
Uniprot ID	TPD52_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 101-151 on human Prostate leucine zipper isoform 5 protein.
Immunogen Region	101-151
Specificity	
Immunogen Sequence	