

TrkA protein (STJ505038)

STJ505038

GENERAL INFORMATION

Product Type	Proteins
Short Description	TrkA-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: IHC: 1:100
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4914
Gene Symbol	NTRK1
Uniprot ID	NTRK1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 360-410 on human TrkA protein.
Immunogen Region	360-410
Specificity	
Immunogen Sequence	