

TRPV6 protein (STJ505041)

STJ505041

GENERAL INFORMATION

Product Type	Proteins
Short Description	TRPV6-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution	
Range	
Formulation	
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthetic peptide taken from the region within amino acids 700-732 of TRPV6 protein. The peptide was covalently modified to achieve desired antigenicity. Modified peptide was covalently coupled to a carrier protein by heterobifunctional cross linker
Immunogen	700-732
Region	
Specificity	
Immunogen	
Sequence	