

C15orf42 protein (STJ505042)

STJ505042

GENERAL INFORMATION

Product Type	Proteins
Short Description	C15orf42-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	90381
Gene Symbol	TICRR
Uniprot ID	TICRR_HUMAN
Immunogen	Synthetic peptide corresponding to unique amino acid sequence on human Treslin protein.
Immunogen Region	
Specificity	
Immunogen Sequence	