

TSH beta protein (STJ505044)

STJ505044

GENERAL INFORMATION

Product Type	Proteins
Short Description	TSH beta-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7252
Gene Symbol	TSHB
Uniprot ID	TSHB_HUMAN
Immunogen	Synthetic peptide corresponding to unique amino acid sequence on human TSH beta chain protein.
Immunogen Region	
Specificity	
Immunogen Sequence	