

TTC27 protein (STJ505046)
STJ505046**GENERAL INFORMATION**

Product Type	Proteins
Short Description	TTC27-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	74196
Gene Symbol	Ttc27
Uniprot ID	TTC27_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 100-150 on tetratricopeptide repeat protein 27 Danio rerio protein.
Immunogen Region	100-150
Specificity	
Immunogen Sequence	