

TYRO3 protein (STJ505048)

STJ505048

GENERAL INFORMATION

Product Type	Proteins
Short Description	TYRO3-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7301
Gene Symbol	TYRO3
Uniprot ID	TYRO3_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 800-850 on human Tyrosine-protein kinase receptor TYRO3 protein.
Immunogen Region	800-850
Specificity	
Immunogen Sequence	