

UFD1L protein (STJ505052)

STJ505052

GENERAL INFORMATION

Product Type	Proteins
Short Description	UFD1L-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7353
Gene Symbol	UFD1
Uniprot ID	UFD1_HUMAN
Immunogen	Purified 20 kDa UFD1L protein
Immunogen Region	
Specificity	
Immunogen Sequence	