

ULK3 protein (STJ505055)

STJ505055

GENERAL INFORMATION

Product Type	Proteins
Short Description	ULK3-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	691171
Gene Symbol	Ulk3
Uniprot ID	ULK3_RAT
Immunogen	Synthetic peptide corresponding to a region within C terminal amino acids 452-472 of Mouse ULK3, coupled to a carrier protein.
Immunogen Region	452-472
Specificity	
Immunogen Sequence	