

VPS26 protein (STJ505062)
STJ505062**GENERAL INFORMATION**

Product Type	Proteins
Short Description	VPS26-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9559
Gene Symbol	VPS26A
Uniprot ID	VP26A_HUMAN
Immunogen	Synthetic peptide taken within amino acids 315-375 on human VPS26 protein.
Immunogen Region	315-375
Specificity	
Immunogen Sequence	