

VPS29 protein (STJ505063)

STJ505063

GENERAL INFORMATION

Product Type	Proteins
Short Description	VPS29-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	51699
Gene Symbol	VPS29
Uniprot ID	VPS29_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 133-183 on rat Vacuolar Protein Sorting-Associated Protein 29.
Immunogen Region	133-183
Specificity	
Immunogen Sequence	