

Visfatin protein (STJ505066)

STJ505066

GENERAL INFORMATION

Product Type	Proteins
Short Description	Visfatin-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: IHC: 1:100
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10135
Gene Symbol	NAMPT
Uniprot ID	NAMPT_HUMAN
Immunogen	C-epitope synthetic peptide taken within amino acid region 400-490 on human Visfatin protein
Immunogen Region	400-490
Specificity	
Immunogen Sequence	