

**WDR5 protein (STJ505068)**  
STJ505068**GENERAL INFORMATION**

<b>Product Type</b>	Proteins
<b>Short Description</b>	WDR5-protein if for use in research applications.
<b>Applications</b>	Immunodepletion, Immunocompetition
<b>Host/Source</b>	
<b>Reactivity</b>	

**PRODUCT PROPERTIES**

<b>Clonality</b>	
<b>Clone ID</b>	
<b>Concentration</b>	
<b>Conjugation</b>	
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	
<b>Isotype</b>	
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

**TARGET INFORMATION**

<b>Gene ID</b>	11091
<b>Gene Symbol</b>	WDR5
<b>Uniprot ID</b>	WDR5_HUMAN
<b>Immunogen</b>	Synthetic peptide: DVKTGKCLKTLPAHSDPV corresponding to amino acids 157-174 of Human WDR5 protein.
<b>Immunogen Region</b>	157-174
<b>Specificity</b>	
<b>Immunogen Sequence</b>	