

MTAP protein (STJ505071)

STJ505071

GENERAL INFORMATION

Product Type	Proteins
Short Description	MTAP-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4507
Gene Symbol	MTAP
Uniprot ID	MTAP_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 233-283 on human methylthioadenosine phosphorylase protein.
Immunogen Region	233-283
Specificity	
Immunogen Sequence	