

XTP4 protein (STJ505079)

STJ505079

GENERAL INFORMATION

Product Type	Proteins
Short Description	XTP4-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	84299
Gene Symbol	MIEN1
Uniprot ID	MIEN1_HUMAN
Immunogen	Multiple epitope synthetic peptide taken within amino acid regions 1-50, 20-70, 40-90 and 80-130 on human XTP4 protein.
Immunogen Region	1-50 20-70 40-90 and 80-130
Specificity	
Immunogen Sequence	