

YES1 protein (STJ505080)

STJ505080

GENERAL INFORMATION

Product Type	Proteins
Short Description	YES1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7525
Gene Symbol	YES1
Uniprot ID	YES_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 450-480 on human YES1 protein.
Immunogen Region	450-480
Specificity	
Immunogen Sequence	