

ZMPSTE24 protein (STJ505082)

STJ505082

GENERAL INFORMATION

Product Type	Proteins
Short Description	ZMPSTE24-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	0.55 µg/µl
Conjugation	
Purification	
Dilution Range	WB: DB: 1:10, 000 ELISA: 1:10, 000
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10269
Gene Symbol	ZMPSTE24
Uniprot ID	FACE1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 215-245 on human ZMPSTE24 protein.
Immunogen Region	215-245
Specificity	
Immunogen Sequence	