

RIP140 protein (STJ505095)
STJ505095**GENERAL INFORMATION**

Product Type	Proteins
Short Description	RIP140-protein if for use in research applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:1, 000
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	268903
Gene Symbol	Nrip1
Uniprot ID	NRIP1_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	