

Soluble Guanylyl Cyclase Alpha 2 protein (STJ505138)

STJ505138

GENERAL INFORMATION

Product Type	Proteins
Short Description	Soluble Guanylyl Cyclase Alpha 2-protein if for use in research applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:750
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	66012
Gene Symbol	Gucy1a2
Uniprot ID	GCYA2_RAT
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	