

## SREBP2 protein (STJ505174)

STJ505174

### GENERAL INFORMATION

|                          |   |
|--------------------------|---|
| <b>Product Type</b>      | Proteins  |
| <b>Short Description</b> | SREBP2-protein if for use in research applications. |
| <b>Applications</b>      | WB  |
| <b>Host/Source</b>       |   |
| <b>Reactivity</b>        |   |

### PRODUCT PROPERTIES

|                            |   |
|----------------------------|---|
| <b>Clonality</b>           |   |
| <b>Clone ID</b>            |   |
| <b>Concentration</b>       |   |
| <b>Conjugation</b>         |   |
| <b>Purification</b>        |   |
| <b>Dilution Range</b>      | WB: 1:500   |
| <b>Formulation</b>         |   |
| <b>Isotype</b>             |   |
| <b>Storage Instruction</b> | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

### TARGET INFORMATION

|                           |             |
|---------------------------|-------------|
| <b>Gene ID</b>            | 6721        |
| <b>Gene Symbol</b>        | SREBF2      |
| <b>Uniprot ID</b>         | SRBP2_HUMAN |
| <b>Immunogen</b>          |             |
| <b>Immunogen Region</b>   |             |
| <b>Specificity</b>        |             |
| <b>Immunogen Sequence</b> |             |