

ETS1 associated protein II protein (STJ505190)

STJ505190

GENERAL INFORMATION

Product Type	Proteins
Short Description	ETS1 associated protein II-protein if for use in research applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:1, 000
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	51567
Gene Symbol	TDP2
Uniprot ID	TYDP2_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	