

TIRAP3b protein (STJ505200)
STJ505200**GENERAL INFORMATION**

Product Type	Proteins
Short Description	TIRAP3b-protein if for use in research applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:250
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	100302736/353376
Gene Symbol	TICAM2
Uniprot ID	TCAM2_HUMAN
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	