

TLR12 protein (STJ505212)

STJ505212

GENERAL INFORMATION

Product Type	Proteins
Short Description	TLR12-protein if for use in research applications.
Applications	WB
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:500
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	384059
Gene Symbol	Tr12
Uniprot ID	TLR12_MOUSE
Immunogen	
Immunogen Region	
Specificity	
Immunogen Sequence	