

**TM6SF2 protein (STJ505215)**  
STJ505215**GENERAL INFORMATION**

|                          |   |
|--------------------------|---|
| <b>Product Type</b>      | Proteins  |
| <b>Short Description</b> | TM6SF2-protein if for use in research applications. |
| <b>Applications</b>      | WB  |
| <b>Host/Source</b>       |   |
| <b>Reactivity</b>        |   |

**PRODUCT PROPERTIES**

|                            |   |
|----------------------------|---|
| <b>Clonality</b>           |   |
| <b>Clone ID</b>            |   |
| <b>Concentration</b>       |   |
| <b>Conjugation</b>         |   |
| <b>Purification</b>        |   |
| <b>Dilution Range</b>      | WB: 1:500   |
| <b>Formulation</b>         |   |
| <b>Isotype</b>             |   |
| <b>Storage Instruction</b> | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

**TARGET INFORMATION**

|                           |                             |
|---------------------------|-----------------------------|
| <b>Gene ID</b>            | <a href="#">107770</a>      |
| <b>Gene Symbol</b>        | <a href="#">Tm6sf2</a>      |
| <b>Uniprot ID</b>         | <a href="#">TM6S2_MOUSE</a> |
| <b>Immunogen</b>          |                             |
| <b>Immunogen Region</b>   |                             |
| <b>Specificity</b>        |                             |
| <b>Immunogen Sequence</b> |                             |