

PDE4D protein (STJ505270)

STJ505270

GENERAL INFORMATION

Product Type	Proteins
Short Description	PDE4D-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5142
Gene Symbol	PDE4B
Uniprot ID	PDE4B_HUMAN
Immunogen	C-terminal cyclic peptide unique to all PDED subtypes variants.
Immunogen Region	C-Term
Specificity	
Immunogen Sequence	