

PDE4B protein (STJ505274)

STJ505274

GENERAL INFORMATION

Product Type	Proteins
Short Description	PDE4B-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5142
Gene Symbol	PDE4B
Uniprot ID	PDE4B_HUMAN
Immunogen	Synthetic peptide common to all PDE4B subtypes taken near N-terminal
Immunogen Region	N-Term
Specificity	
Immunogen Sequence	