

## PDE4B protein (STJ505275)

STJ505275

### GENERAL INFORMATION

<b>Product Type</b>	Proteins
<b>Short Description</b>	PDE4B-protein if for use in research applications.
<b>Applications</b>	WB
<b>Host/Source</b>	
<b>Reactivity</b>	

### PRODUCT PROPERTIES

<b>Clonality</b>	
<b>Clone ID</b>	
<b>Concentration</b>	
<b>Conjugation</b>	
<b>Purification</b>	
<b>Dilution Range</b>	WB: 1:500-1:1, 000
<b>Formulation</b>	
<b>Isotype</b>	
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	5142
<b>Gene Symbol</b>	PDE4B
<b>Uniprot ID</b>	PDE4B_HUMAN
<b>Immunogen</b>	
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	