

TICAM2 protein (STJ505282)

STJ505282

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Proteins |
| Short Description | TICAM2-protein if for use in research applications. |
| Applications | Immunodepletion, Immunocompetition |
| Host/Source | |
| Reactivity | |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | |
| Clone ID | |
| Concentration | |
| Conjugation | |
| Purification | |
| Dilution Range | |
| Formulation | |
| Isotype | |
| Storage Instruction | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|--|
| Gene ID | 100302736/353376 |
| Gene Symbol | TICAM2 |
| Uniprot ID | TCAM2_HUMAN |
| Immunogen | Synthetic peptide corresponding to unique amino acid sequences on TRIF-related adapter molecule. |
| Immunogen Region | |
| Specificity | |
| Immunogen Sequence | |