

WVOX 1 + 4 protein (STJ505284)

STJ505284

GENERAL INFORMATION

Product Type	Proteins
Short Description	WVOX 1 + 4-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	51741
Gene Symbol	WVOX
Uniprot ID	WVOX_HUMAN
Immunogen	Synthetic peptide corresponding to unique amino acid sequence on WVOX 1 and WVOX4 protein.
Immunogen Region	
Specificity	
Immunogen Sequence	