

TOMM40 protein (STJ505296)

STJ505296

GENERAL INFORMATION

Product Type	Proteins
Short Description	TOMM40-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10452
Gene Symbol	TOMM40
Uniprot ID	TOM40_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 300-350 on human Mitochondrial import receptor subunit TOM40 homolog.
Immunogen Region	300-350
Specificity	
Immunogen	
Sequence	