

DCK protein (STJ505299)

STJ505299

GENERAL INFORMATION

Product Type	Proteins
Short Description	DCK-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	13178
Gene Symbol	Dck
Uniprot ID	DCK_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 210-260 on mouse Deoxycytidine kinase.
Immunogen Region	210-260
Specificity	
Immunogen Sequence	