

Ephrin Receptor A1 protein (STJ505303)

STJ505303

GENERAL INFORMATION

Product Type	Proteins
Short Description	Ephrin Receptor A1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2041
Gene Symbol	EPHA1
Uniprot ID	EPHA1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 55-605 on human Ephrin type-A receptor 1.
Immunogen Region	55-605
Specificity	
Immunogen Sequence	