

Insulin Receptor protein (STJ505308)

STJ505308

GENERAL INFORMATION

Product Type	Proteins
Short Description	Insulin Receptor-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	WB: 1:500
Formulation	
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	3643
Gene Symbol	INSR
Uniprot ID	INSR_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 900-950 on human Insulin receptor, transcript variant 2 protein. Cross reactive with both Insulin receptor alpha and insulin receptor beta variants.
Immunogen Region	900-950
Specificity	
Immunogen Sequence	