

Insulin Receptor beta protein (STJ505311)

STJ505311

GENERAL INFORMATION

Product Type	Proteins
Short Description	Insulin Receptor beta-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3643
Gene Symbol	INSR
Uniprot ID	INSR_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 725-775 on human Insulin receptor beta subunit protein.
Immunogen Region	725-775
Specificity	
Immunogen Sequence	