

MAP3K4 protein (STJ505315)

STJ505315

GENERAL INFORMATION

Product Type	Proteins
Short Description	MAP3K4-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4216
Gene Symbol	MAP3K4
Uniprot ID	M3K4_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1558-1608 on human MAP3K4 protein.
Immunogen Region	1558-1608
Specificity	
Immunogen Sequence	