

ALDH1A1 protein (STJ505317)

STJ505317

GENERAL INFORMATION

Product Type	Proteins
Short Description	ALDH1A1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	24188
Gene Symbol	Aldh1a1
Uniprot ID	AL1A1_RAT
Immunogen	Synthetic peptide taken within amino acid region 441-491 on rat retinal dehydrogenase 1 protein.
Immunogen Region	441-491
Specificity	
Immunogen Sequence	