

**ALDH1A1 protein (STJ505318)**  
STJ505318**GENERAL INFORMATION**

|                          |  |
|--------------------------|--|
| <b>Product Type</b>      | Proteins   |
| <b>Short Description</b> | ALDH1A1-protein if for use in research applications. |
| <b>Applications</b>      | WB   |
| <b>Host/Source</b>       |  |
| <b>Reactivity</b>        |  |

**PRODUCT PROPERTIES**

|                            |   |
|----------------------------|---|
| <b>Clonality</b>           |   |
| <b>Clone ID</b>            |   |
| <b>Concentration</b>       |   |
| <b>Conjugation</b>         |   |
| <b>Purification</b>        |   |
| <b>Dilution Range</b>      | WB: 1:2, 500-1:5, 000   |
| <b>Formulation</b>         |   |
| <b>Isotype</b>             |   |
| <b>Storage Instruction</b> | Store at -20°C for long term storage. Avoid freeze-thaw cycles. |

**TARGET INFORMATION**

|                           |                           |
|---------------------------|---------------------------|
| <b>Gene ID</b>            | <a href="#">24188</a>     |
| <b>Gene Symbol</b>        | <a href="#">Aldh1a1</a>   |
| <b>Uniprot ID</b>         | <a href="#">AL1A1_RAT</a> |
| <b>Immunogen</b>          |                           |
| <b>Immunogen Region</b>   |                           |
| <b>Specificity</b>        |                           |
| <b>Immunogen Sequence</b> |                           |