

VEGFR1 protein (STJ505319)

STJ505319

GENERAL INFORMATION

Product Type	Proteins
Short Description	VEGFR1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2321
Gene Symbol	FLT1
Uniprot ID	VGFR1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 1250-1300 on human Vascular endothelial growth factor receptor 1 protein.
Immunogen Region	1250-1300
Specificity	
Immunogen	
Sequence	