

## Estrogen Receptor alpha protein (STJ505322)

STJ505322

### GENERAL INFORMATION

<b>Product Type</b>	Proteins
<b>Short Description</b>	Estrogen Receptor alpha-protein if for use in research applications.
<b>Applications</b>	Immunodepletion, Immunocompetition
<b>Host/Source</b>	
<b>Reactivity</b>	

### PRODUCT PROPERTIES

<b>Clonality</b>	
<b>Clone ID</b>	
<b>Concentration</b>	
<b>Conjugation</b>	
<b>Purification</b>	
<b>Dilution Range</b>	
<b>Formulation</b>	
<b>Isotype</b>	
<b>Storage Instruction</b>	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	2099
<b>Gene Symbol</b>	ESR1
<b>Uniprot ID</b>	ESR1_HUMAN
<b>Immunogen</b>	Synthetic peptide taken within amino acid region 520-570 on human Estrogen receptor alpha protein.
<b>Immunogen Region</b>	520-570
<b>Specificity</b>	
<b>Immunogen Sequence</b>	