

Sialoadhesin protein (STJ505324)

STJ505324

GENERAL INFORMATION

Product Type	Proteins
Short Description	Sialoadhesin-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	0.5 µg/µl
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6614
Gene Symbol	SIGLEC1
Uniprot ID	SN_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 10-60 on human Sialoadhesin protein.
Immunogen Region	10-60
Specificity	
Immunogen Sequence	