

Humanin protein (STJ505326)

STJ505326

GENERAL INFORMATION

Product Type	Proteins
Short Description	Humanin-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	NA
Gene Symbol	MT-RNR2
Uniprot ID	HUNIN_HUMAN
Immunogen	Synthetic peptide corresponding to amino acid region 5-24 on Humanin protein. Sequence: GFSCLLSTTATDLPVKRRT
Immunogen Region	
Specificity	
Immunogen Sequence	