

ENT1 protein (STJ505328)

STJ505328

GENERAL INFORMATION

Product Type	Proteins
Short Description	ENT1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage	Store at -20°C for long term storage. Avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	63959
Gene Symbol	Slc29a1
Uniprot ID	S29A1_MOUSE
Immunogen	Synthetic peptide taken within amino acid region 410-460 on mouse Equilibrative nitrobenzylthioinosine-sensitive nucleoside transporter ENT1 protein.
Immunogen Region	410-460
Specificity	
Immunogen Sequence	