

SERINC1 protein (STJ505330)

STJ505330

GENERAL INFORMATION

Product Type	Proteins
Short Description	SERINC1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	57515
Gene Symbol	SERINC1
Uniprot ID	SERC1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 250-310 on human Serine incorporator 1 protein.
Immunogen Region	250-310
Specificity	
Immunogen Sequence	