

FPR1 protein (STJ505336)
STJ505336**GENERAL INFORMATION**

Product Type	Proteins
Short Description	FPR1-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2357
Gene Symbol	FPR1
Uniprot ID	FPR1_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 150-200 on N-formyl peptide receptor.
Immunogen Region	150-200
Specificity	
Immunogen Sequence	