

PDE2A2 protein (STJ505340)

STJ505340

GENERAL INFORMATION

Product Type	Proteins
Short Description	PDE2A2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5138
Gene Symbol	PDE2A
Uniprot ID	PDE2A_HUMAN
Immunogen	Synthetic peptide corresponding to region within amino acids 1-50 on N-terminus of human PDE2A2 protein.
Immunogen Region	1-50aa
Specificity	
Immunogen Sequence	