

Pax2 protein (STJ505344)

STJ505344

GENERAL INFORMATION

Product Type	Proteins
Short Description	Pax2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5076
Gene Symbol	PAX2
Uniprot ID	PAX2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 255-305 on human PAX2 protein.
Immunogen Region	255-305
Specificity	
Immunogen Sequence	