

Inhibin alpha protein (STJ505346)

STJ505346

GENERAL INFORMATION

Product Type	Proteins
Short Description	Inhibin alpha-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3623
Gene Symbol	INHA
Uniprot ID	INHA_HUMAN
Immunogen	Synthetic peptides corresponding to amino acid region taken within 1-50 and 220-270 on Inhibin alpha gene.
Immunogen Region	1-50 and 220-270
Specificity	
Immunogen Sequence	