

SPACA3 protein (STJ505347)

STJ505347

GENERAL INFORMATION

Product Type	Proteins
Short Description	SPACA3-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	124912
Gene Symbol	SPACA3
Uniprot ID	SACA3_HUMAN
Immunogen	Synthetic peptides taken within amino acid region 20-70 on human sperm acrosome membrane-associated protein 3 isoform 2.
Immunogen Region	20-70
Specificity	
Immunogen	
Sequence	