

STAT3 protein (STJ505353)
STJ505353**GENERAL INFORMATION**

Product Type	Proteins
Short Description	STAT3-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6774
Gene Symbol	STAT3
Uniprot ID	STAT3_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 700-750 on human STAT3 protein.
Immunogen Region	700-750
Specificity	
Immunogen Sequence	