

FGFR2 protein (Phospho) (STJ505361)

STJ505361

GENERAL INFORMATION

Product Type	Proteins
Short Description	Phospho-FGFR2-protein if for use in research applications.
Applications	Immunodepletion, Immunocompetition
Host/Source	
Reactivity	

PRODUCT PROPERTIES

Clonality	
Clone ID	
Concentration	
Conjugation	
Purification	
Dilution Range	
Formulation	
Isotype	
Storage Instruction	Store at -20°C for long term storage. Avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2263
Gene Symbol	FGFR2
Uniprot ID	FGFR2_HUMAN
Immunogen	Synthetic peptide taken within amino acid region 520-570 on human Fibroblast Growth Factor receptor 2 isoform X6 protein.
Immunogen Region	520-570
Specificity	
Immunogen	
Sequence	